


# The programme button

- Your hearing aid can have up to 4 different programmes. The programmes can be selected using the button at the top, on the back of your hearing aid.
- To change programmes press the button. The number of beeps you hear indicates the programme you are in.  
1 Beep = Programme 1  
2 Beeps = Programme 2, etc.
- The programmes cycle through 1 to 2 to 3 to 4, then back to 1.
- When you switch the hearing aid off and back on again it will always return to programme 1.



## Record your programmes here

	Programme Name	When to use
<b>Programme 1</b>		
<b>Programme 2</b>		
<b>Programme 3</b>		
<b>Programme 4</b>		

Your hearing aid has been set up to beep when you change between programmes.

Yes     No

Audiologist contact details:

---

The right choice for the NHS

ReSound Up™

UK 05.63 A 7.2017